

Single Use Channel Cleaning Brush


BW-201B/T

Single Use Channel-Opening Cleaning Brush

MAJ-1339



Single Use Olympus Cleaning Brushes



As an endoscope manufacturer, Olympus is pleased to offer a brush that has been designed to facilitate maximum cleaning performance while ensuring compatibility with the endoscope channel.

Specifications

	BW-201B	BW-201T
Compatible Channel Diameter	2~3.2mm	2~4.2 mm
Working Length	950 mm	2200 mm

BW-201B/T.A set

Single Use Channel Cleaning Brush

Designed for one-time use, this brush is ideal for cleaning out suction channels in compatible Olympus endoscopes and openings in compatible Endo-Therapy accessories. Individually packed — 50 packs per box.



BW-201B/T.B set

Channel Cleaning Brush and Channel-Opening Cleaning Brush


This set contains both the Single Use channel cleaning brush and the Single Use channel-opening cleaning brush. One of each per pack — 50 packs per box.



MAJ-1339

Single Use Channel-Opening Cleaning Brush

Specifically designed for cleaning scope channel ports. Single Use only. Individually packed — 50 packs per box.



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer



OLYMPUS MEDICAL SYSTEMS CORP. OLYMPUS MEDICAL SYSTEMS EUROPA GMBH OLYMPUS AMERICA INC. OLYMPUS LATIN AMERICA, INC.

a 33126-2097, U.S.A 5301 Blue Lagoon Drive, Sont Table KEYMED LTD.
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD. OLYMPUS HONG KONG AND CHINA LIMITED OLYMPUS (BEIJING) SALES & SERVICE CO.,LTD. OLYMPUS MOSCOW LIMITED LIABILITY COMPANY OLYMPUS AUSTRALIA PTY LTD