


**OLYMPUS**<sup>®</sup>

Your Vision, Our Future

OLYMPUS  
**EndoTherapy**  
BEYOND THE SCOPE

● EndoTherapy for small intestinal procedures



# Versatile EndoTherapy accessories from Olympus offer exceptional support for small intestinal procedures, advancing enteroscopy into the future


## Sampling

### Biopsy Forceps



#### Single Use

Employment of an improved sheath section enables smooth insertion into a scope channel. These forceps are especially designed to prevent damage to the interior of a scope channel. Moreover, the newly developed blades (double-step blades) on the cups ensure sharp cutting quality, and swing jaw mechanism enables easier cup positioning for biopsy sites approached tangentially.



- FB-210U** Alligator Jaw without needle
- FB-220U** Alligator Jaw with needle
- FB-230U** Standard oval without needle
- FB-240U** Standard oval with needle

Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	20 pcs / Box

#### Reusable

A hole of approximately 1 mm diameter at the center of the cups minimizes tissue damage and provides large sample yield.



- FB-28Y-1** Standard oval without needle
- FB-28Z-1** Standard oval without needle


Minimum Channel Size	2.8 mm
Working Length	FB-28Y-1 : 2700 mm FB-28Z-1 : 3000 mm
Quantity	1 pc / Box

## Foreign Body Retrieval

### Retrieval Basket

Especially designed for foreign body removal in small intestine, this retrieval basket offers reliable performance in retrieving resected tissue and other foreign objects with 6 wire basket and rotatable function.

#### Single Use



- FG-460YR** 6 wire basket

Minimum Channel Size	2.8 mm
Working Length	2700 mm
Basket Design	6 wire basket
Rotatability	Yes
Quantity	1 pc / Box


## Polypectomy

### Hot Biopsy Forceps



Disposable hot biopsy forceps with the slim outer diameter, the proprietary sheath design and improved flexibility ensure smooth insertion into the endoscope channel.

#### Single Use



- FD-210U** Alligator Jaw
- FD-230U** Standard oval


Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	5 pcs / Box

### Electrosurgical Snare



#### Single Use

Disposable design with integrated handle for ease of use.




- SD-221U-25** Crescent type

Minimum Channel Size	2.0 mm
Working Length	2300 mm
Loop Opening Width	25 mm
Quantity	10 pcs / Box

#### Reusable

Especially designed for small intestinal procedure.



- SD-310Z-25** Oval type

Minimum Channel Size	2.8 mm
Working Length	3000 mm
Loop Opening Width	25 mm
Quantity	2 Wires and 2 tube sheaths

With the introduction of an innovative, yet surprisingly simple Single Balloon Enteroscope System, Olympus has literally opened up the small intestine to enteroscopy, helping this promising field take a quantum leap forward. Our sophisticated line of versatile EndoTherapy accessories combine with the advanced capabilities of the Single Balloon Enteroscope to create a synergy that will put you at the cutting edge of enteroscopy.

## Hemostasis

### Injection Needle

Clear catheter provides proper flexibility and rigidity. Especially suitable for submucosal injection in hemostasis and EMR.



#### Single Use



**NM-200U-0423** Short bevel

Minimum Channel Size	2.8 mm
Working Length	2300 mm
Needle Length	4 mm
Gauge	23 G
Quantity	6 pcs / Box

## Other Accessories

### Washing Pipe / Spray Catheters

Especially designed for spraying dye, these catheters enable circular spraying of the liquid from the distal end.

#### Reusable



**PW-1V-1** Standard type

**PW-5V-1** Spray type


Minimum Channel Size	2.8 mm
Working Length	2400 mm
Feature	PW-1V-1 : Inside Chip (Routine washing) PW-5V-1 : Screw tip (Staining)
Quantity	PW-1V-1 : 2 pcs / Box (A), 1 pc / Box (B) PW-5V-1 : 1 pc / Box

### Cleaning Brushes



Designed for one-time use, this brush is ideal for cleaning suction channels and openings of compatible Olympus endoscopes.

#### Single Use




**BW-201T(A)** Single use channel cleaning brush

**BW-201T(B)** Channel cleaning brush and channel-opening cleaning brush

Working Length	2200 mm
Quantity	A : 50 pcs / Box B : 50 sets / Box

#### Reusable



**BW-9Y** Coil sheath channel cleaning brush


**MH-507** Channel-port cleaning brush

Working Length	2850 mm
Quantity	1 pc / box

### Distal Attachment

When attached to the distal end of a scope, this transparent cap assists in observation of the mucosa and helps achieve a good field of view.

#### Single Use



**D-201-10704** Soft Cap, Straight type

Outer Diameter	10.7 mm
Distal Length from Scope	4 mm
Applicable Scope	SIF-Q180, Q260
Quantity	12 pcs / box



OLYMPUS

Endo**Therapy**  
BEYOND THE SCOPE

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**<sup>®</sup>

**OLYMPUS MEDICAL SYSTEMS CORP.**  
Shinjuku Monolith, 3-1 Nishi-Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914, Japan  
**OLYMPUS MEDICAL SYSTEMS EUROPA GMBH**  
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
3500 Corporate Parkway, P.O. Box 610, Center Valley, PA18034-0610, U.S.A.  
**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, Florida 33126-2097, U.S.A.  
**KEYMED LTD.**  
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373  
**OLYMPUS HONG KONG AND CHINA LIMITED**  
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong  
**OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.**  
Room.1202, NCI Tower, A12 Jianguomenwai Avenue Chaoyang District Beijing, 100022, China  
**OLYMPUS MOSCOW LIMITED LIABILITY COMPANY**  
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia  
**OLYMPUS AUSTRALIA PTY LTD**  
31 Gilby Road, Mount Waverley, VIC., 3149, Australia

[www.olympus.com](http://www.olympus.com)

F1228TT-0707