


OLYMPUS[®]

Your Vision, Our Future

OLYMPUS
EndoTherapy
BEYOND THE SCOPE

● EndoTherapy for small intestinal procedures



Versatile EndoTherapy accessories from Olympus offer exceptional support for small intestinal procedures, advancing enteroscopy into the future

Sampling

Biopsy Forceps



Single Use

Employment of an improved sheath section enables smooth insertion into a scope channel. These forceps are especially designed to prevent damage to the interior of a scope channel. Moreover, the newly developed blades (double-step blades) on the cups ensure sharp cutting quality, and swing jaw mechanism enables easier cup positioning for biopsy sites approached tangentially.



- FB-210U** Alligator Jaw without needle
- FB-220U** Alligator Jaw with needle
- FB-230U** Standard oval without needle
- FB-240U** Standard oval with needle

Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	20 pcs / Box

Reusable

A hole of approximately 1 mm diameter at the center of the cups minimizes tissue damage and provides large sample yield.



- FB-28Y-1** Standard oval without needle
- FB-28Z-1** Standard oval without needle


Minimum Channel Size	2.8 mm
Working Length	FB-28Y-1 : 2700 mm FB-28Z-1 : 3000 mm
Quantity	1 pc / Box

Foreign Body Retrieval

Retrieval Basket

Especially designed for foreign body removal in small intestine, this retrieval basket offers reliable performance in retrieving resected tissue and other foreign objects with 6 wire basket and rotatable function.

Single Use



- FG-460YR** 6 wire basket

Minimum Channel Size	2.8 mm
Working Length	2700 mm
Basket Design	6 wire basket
Rotatability	Yes
Quantity	1 pc / Box


Polypectomy

Hot Biopsy Forceps



Disposable hot biopsy forceps with the slim outer diameter, the proprietary sheath design and improved flexibility ensure smooth insertion into the endoscope channel.

Single Use



- FD-210U** Alligator Jaw
- FD-230U** Standard oval


Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	5 pcs / Box

Electrosurgical Snare



Single Use

Disposable design with integrated handle for ease of use.




- SD-221U-25** Crescent type

Minimum Channel Size	2.0 mm
Working Length	2300 mm
Loop Opening Width	25 mm
Quantity	10 pcs / Box

Reusable

Especially designed for small intestinal procedure.



- SD-310Z-25** Oval type

Minimum Channel Size	2.8 mm
Working Length	3000 mm
Loop Opening Width	25 mm
Quantity	2 Wires and 2 tube sheaths

With the introduction of an innovative, yet surprisingly simple Single Balloon Enteroscope System, Olympus has literally opened up the small intestine to enteroscopy, helping this promising field take a quantum leap forward. Our sophisticated line of versatile EndoTherapy accessories combine with the advanced capabilities of the Single Balloon Enteroscope to create a synergy that will put you at the cutting edge of enteroscopy.

Hemostasis

Injection Needle

Clear catheter provides proper flexibility and rigidity. Especially suitable for submucosal injection in hemostasis and EMR.



Single Use



NM-200U-0423 Short bevel

Minimum Channel Size	2.8 mm
Working Length	2300 mm
Needle Length	4 mm
Gauge	23 G
Quantity	6 pcs / Box

Other Accessories

Washing Pipe / Spray Catheters

Especially designed for spraying dye, these catheters enable circular spraying of the liquid from the distal end.

Reusable



PW-1V-1 Standard type

PW-5V-1 Spray type


Minimum Channel Size	2.8 mm
Working Length	2400 mm
Feature	PW-1V-1 : Inside Chip (Routine washing) PW-5V-1 : Screw tip (Staining)
Quantity	PW-1V-1 : 2 pcs / Box (A), 1 pc / Box (B) PW-5V-1 : 1 pc / Box

Cleaning Brushes



Designed for one-time use, this brush is ideal for cleaning suction channels and openings of compatible Olympus endoscopes.

Single Use




BW-201T(A) Single use channel cleaning brush

BW-201T(B) Channel cleaning brush and channel-opening cleaning brush

Working Length	2200 mm
Quantity	A : 50 pcs / Box B : 50 sets / Box

Reusable



BW-9Y Coil sheath channel cleaning brush


MH-507 Channel-port cleaning brush

Working Length	2850 mm
Quantity	1 pc / box

Distal Attachment

When attached to the distal end of a scope, this transparent cap assists in observation of the mucosa and helps achieve a good field of view.

Single Use



D-201-10704 Soft Cap, Straight type

Outer Diameter	10.7 mm
Distal Length from Scope	4 mm
Applicable Scope	SIF-Q180, Q260
Quantity	12 pcs / box



OLYMPUS

Endo**Therapy**
BEYOND THE SCOPE

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS[®]

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 3-1 Nishi-Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914, Japan
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
3500 Corporate Parkway, P.O. Box 610, Center Valley, PA18034-0610, U.S.A.
OLYMPUS LATIN AMERICA, INC.
5301 Blue Lagoon Drive, Suite 290 Miami, Florida 33126-2097, U.S.A.
KEYMED LTD.
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373
OLYMPUS HONG KONG AND CHINA LIMITED
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.
Room.1202, NCI Tower, A12 Jianguomenwai Avenue Chaoyang District Beijing, 100022, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
OLYMPUS AUSTRALIA PTY LTD
31 Gilby Road, Mount Waverley, VIC., 3149, Australia

www.olympus.com

F1228TT-0707